

Protein A/G Magnetic Beads

CATALOG #: 6527-1

AMOUNT: 1 ml

LOT #: _____

PREPARATION: Protein A/G Magnetic Beads are prepared by covalently coupling Recombinant fusion Protein A/G (contains eight IgG binding domains, BV catalog # 6502) to 6% crosslinked magnetically beaded agarose. The coupling technique is optimized to give a high binding capacity for IgG. The capacity of IgG binding is generally greater than 10 mg of human IgG per ml of wet gel.

CONTENTS: Supplied as a 50% slurry in PBS with 0.02% sodium azide.

TECHNICAL SPECIFICATIONS:

Parameter	Description
Support Characteristics	Paramagnetic, spherical, 6 % cross-linked agarose
Ligand	Recombinant fusion Protein A/G
Particle Size	75 – 150 μ m
Binding Capacity	Generally >10 mg human IgG/ml wet beads
Working Temperature	Room temperature
Storage Solution	PBS w/0.02% NaN ₃
Storage Temperature	4 – 8 °C
Stability	Stable, as supplied, for at least 1 year.

FEATURES: Easy to use, high-binding capacity, non-adherent beads. Useful for immunoprecipitation and enrichment of IgG antibodies. High affinity for Fc region of IgG antibodies from a variety of species. Protein A/G binds to all IgG subclasses from various mammalian species, including all IgGs that bind to both Protein A and Protein G individually, making it the ideal choice for purification of all kinds of polyclonal or monoclonal IgG antibodies.

USAGE: For Research Purpose Only! Not to be used in humans!

Suggested Protocol:

Prepare the antibody solution by diluting the required amount of antibody in binding buffer before running the protocol.

- Magnetic Bead Preparation (perform three times)
 - Dispense the required amount of magnetic beads into a 1.5 ml microfuge tube.
 - Place the tube in the magnetic rack and remove the storage solution.
 - Add 500 μ l binding buffer.
 - Resuspend the beads.
 - Remove the liquid
- Antibody Capture
 - Immediately add the antibody solution.
 - Resuspend and mix (slow end-over-end) for at least 15 minutes.
 - Remove the liquid.
- Washing
 - Add 500 μ l Binding Buffer containing 0.5 M NaCl; Remove the liquid.
 - Add 500 μ l Binding Buffer; Remove the liquid.
- Target Binding
 - Add sample diluted in binding buffer.
 - Incubate with slow end-over-end mixing for up to 60 minutes.
 - Remove and collect unbound fraction.
- Washing (perform three times)
 - Add 500 μ l wash buffer
 - Remove liquid (save washes to troubleshoot)
- Elution (perform three times)
 - Add 2 volumes elution buffer (vs. bead volume).
 - Completely resuspend beads and incubate at least 2 minutes.
 - Remove and collect elution fraction.

Recommended Buffer Examples:

Binding buffer: 50 mM Tris, 150 mM NaCl, pH 7.5
Wash buffer: 50 mM Tris, 150 mM NaCl, pH 7.5 (or add 1% Octylglucoside to this buffer)
 (Could also try 1X PBS as both binding and wash buffer)
Elution buffer: 0.1 M -0.2 M Glycine pH 2.5-3.1 (or 0.1 M citric acid, pH 2.5-3.1 or 2.5 % Acetic Acid)

RELATED PRODUCTS:

Recombinant Protein A	Protein A Sepharose	Protein A Magnetic Beads
Recombinant Protein G	Protein G Sepharose	Protein G Magnetic Beads
Recombinant Protein L	Protein L Sepharose	Protein L Magnetic Beads
Recombinant Protein A/G	Protein A/G Sepharose	Protein A/G Magnetic Beads
Recombinant Protein A/G/L	Protein A/G/L Sepharose	Protein A/G/L Magnetic Beads
Protein G Polyclonal Antibody	Protein G-Biotin	Protein G-FITC
Protein A Polyclonal Antibody	Protein G Coated Plate	
Protein L Polyclonal Antibody		