

## Laq™ DNA Polymerase

(2 Units/μl; Store at -20°C)

<b>CATALOG #:</b>	9004-500	500 units
	9004-2500	2500 units

### DESCRIPTION:

Thermostable DNA polymerase of *Pyrococcus glycovorans* from the East Pacific Rise has been isolated and characterized. It amplifies **L**ong DNA products with high **A**ccurate and **Q**uick polymerization rate in DNA Polymerization Chain Reaction (PCR). We named the DNA Polymerase Laq™ DNA Polymerase. The Laq™ DNA polymerase provides the most accurate, consistent performance with high-fidelity up to mutation frequency as low as 10<sup>-7</sup>-10<sup>-6</sup>. The high processivity of the Laq™ DNA Polymerase shortens PCR cycle time to 1/3-1/2 of regular polymerases and amplify up to 10 or 12kb DNA products. It can be used in most protocol that requires high fidelity, quick amplification and long PCR Reactions.

### COMPONENTS:

Component	9004-500	9004-2500
Laq™ DNA Polymerase	500 units	2500 units
Opti-PCR Buffer I	1 x 2 ml	5 x 2 ml

**CONCENTRATION:** 2 Units/μl

### UNIT DEFINITION:

One unit is defined as the amount of enzyme required to catalyze the incorporation of 10 nanomoles of dNTPs into acid insoluble material in 30 minutes at 74°C under standard DNA polymerase assay conditions.

### ENZYME STORAGE BUFFER:

20 mM Tris-HCl, pH 8.0, 100 mM KCl, 0.1 mM EDTA, 1 mM DTT, 50% glycerol, 0.5% Tween 20 and 200 μg/ml BSA.

### 10X PCR BUFFER:

- **Opti-PCR Buffer I:** Used for high fidelity long PCR cloning.
- **Opti-PCR Buffer II:** Used for long PCR products and/or high GC-rich template.
- The Laq™ polymerase also works with most of other polymerase reaction buffers.

### SUGGESTED PCR REACTION MIX:

For each 25 μl PCR reaction, mix the following components in a thin-walled PCR tube:


- |                                      |               |
|--------------------------------------|---------------|
| • 10X PCR Buffer                     | 2.5 μl        |
| • dNTP Mix (2.5 mM each dNTP)        | 2.0 μl        |
| • Laq™ DNA Polymerase (2.0 u/μl)     | 0.2-1.0 units |
| • Template DNA                       | 0.1-500 ng    |
| • Upstream Primer                    | 50-100 pmol   |
| • Downstream Primer                  | 50-100 pmol   |
| • H <sub>2</sub> O to a total volume | 25 μl         |

### SUGGESTED PCR CYCLES:

Initial Denaturation:	98°C	30 sec.
Denaturation:	98°C	10 sec.
Annealing:	50-72°C	5-10 sec.
Elongation:	72°C	15-30 sec/kb.
Final Elongation:	72°C	7 min.

### Note:

For optimal specificity of GC-rich template DNA or expected PCR product longer than 10kb, further optimizing the concentration of MgCl<sub>2</sub>, temperature, cycling time and adding appropriate DMSO could improve PCR performance.



**FIGURE LEGEND:** DNA polymerase comparisons: PCR reactions are performed using Taq, Pfu, Phusion and Laq™ DNA polymerase in their optimized reaction buffers. Laq™ DNA Polymerase can routinely amplify 8 kb DNA fragment, and up to 12 kb DNA fragment when optimizing the reaction conditions.

### RELATED PRODUCTS:

- InsertFinder Insert Quick Screening Kit
- Mitochondrial DNA Isolation Kit
- Luciferase Reporter Assay Kit
- T4 DNA Ligase
- dNTP Mix
- Gel-FAST™ 20 Minutes Gel Staining/Destaining Kit
- Protease Inhibitor cocktail
- Link-FAST™ 5 Minutes DNA Ligation Kit
- Genomic DNA Isolation Kit
- β-Galactosidase Staining Kit
- Agarase
- siRNA Vectors

**FOR RESEARCH USE ONLY! Not to be used in humans**