

ToxOut™ Protein A (Sephacrose) Column

(Catalog# M2503- 1, 5; 1, 5 columns, Store at 4°C)

11/18

I. Introduction:

Protein A beads have been widely used for IgG purification. BioVision's Protein A (Cat# 6500) is a genetically engineered protein containing five IgG-binding regions of native Protein A. Maximum IgG specificity is achieved by eliminating its cell wall binding region, albumin binding region, and other non-specific regions. ToxOut™ Endotoxin-Free Protein A Sepharose beads (Cat#M1300) for IgG purification contain covalently coupled recombinant Protein A to 6% cross-linked Sepharose beads, the most popular resin for protein purification. ToxOut™ Endotoxin-Free Protein A (Sephacrose) column is designed for endotoxin-free IgG purification and immunoprecipitation. This column can rapidly purify IgG from serum, ascites and cell culture media from various species such as human, mouse, rat, goat and rabbit.

II. Application:

- Purify up to 20 mg/ml IgG from endotoxin-free serum, ascites and cell culture media.

III. Contents:

M2503-1 ToxOut™ Protein A (Sephacrose) Column Qty: 1

M2503-5 ToxOut™ Protein A (Sephacrose) Column Qty: 5

IV. Storage and Handling:



- Store at 4°C. Columns expire 2 months after opening.
- Read the entire protocol before performing the assay.
- To prevent endotoxin contamination from dust, solution or dirty lab ware, only use endotoxin-free solutions and tubes and proceed with extra caution.
- Do not let the resin dry out anytime.

V. Specification:

ToxOut™ Protein A (Sephacrose) Column capacity: ~20mg/ml

| Species | IgG Binding Strength |
|---------|----------------------|
| Human | ++++ |
| Mouse | ++++ |
| Rabbit | ++++ |
| Rat | + |
| Goat | + |

Note: Strong Binding +++++/ Weak Binding+



Figures a) Endotoxin level in concentration (left) and total amount (right) before and after rabbit IgG were purified by ToxOut™ Endotoxin Free Protein A (Sephacrose) Column. The recovered IgG shows > 96% endotoxin reduction. **b)** Elution curve for rabbit IgG purification using this Column. The IgG recovery rate from first three elution fraction is 91.5%.

FOR RESEARCH USE ONLY! Not to be used on humans.

FOR RESEARCH USE ONLY!



Related Products:

- ToxOut™ Rapid Endotoxin Removal Kit (K2501)
- ToxOut™ High Capacity Endotoxin Removal Kit (K2502)
- ToxOut™ Protein A-Sepharose Antibody Purification Kit (K2503)
- ToxOut™ Endotoxin Free Protein A (Sepharose) beads (M1300)
- ToxOut™ Endotoxin Free Protein G (Sepharose) beads (M1301)
- ToxOut™ Endotoxin Free Collection Tube (7937)
- Protein A (Cat# 6500, 6520)

FOR RESEARCH USE ONLY! Not to be used on humans.

155 S. Milpitas Blvd., Milpitas, CA 95035 USA | T: (408)493-1800 F: (408)493-1801 | www.biovision.com | tech@biovision.com