## **BioVision**

## PBS Tablets (1000 ml)

**ALTERNATE NAME:** Phosphate Buffered Saline Tablets

**CATALOG #**: 2129-10, 100

**AMOUNT:** 10 Tablets, 100 Tablets

APPEARANCE: Solid

FORMULATION: Dissolving one tablet in 1 liter of deionized H<sub>2</sub>O yields a

solution containing 140 mM NaCl, 10 mM phosphate buffer,

and 2.7 mM KCl, pH 7.4 at 25°C.

**RECONSTITUTION:** Dissolve one tablet in one liter deionized H₂O to yield 10 mM

phosphate buffer, pH 7.4, 140 mM NaCl, 2.7 mM KCl.

**STORAGE:** Store at room temperature. After reconstitution, store at +4 °C

**DESCRIPTION:** One of the most commonly used buffers provided in a

convenient, tablet form. PBS is isotonic and non-toxic to cells and has the ability to maintain their osmolarity thereby making the buffer suitable for washing procedures in cell cultures and

for immunoassays.

rev.09/15

RELATED PRODUCTS:

ASB-14 (Cat. No. 1738-5G, 25G)

ASB-16 (Cat. No. 2136-1G, 5G)

CHAPS, High Purity (Cat. No. 1545-10G, 50G, 500G, 1 Kg) Digitonin, Water-soluble (Cat. No. 2081-100, 250, 1000)

EZSolution™ Digitonin, Water-soluble (Cat. No. 2082-1)

NP-40 Substitute Detergent, MegaPure™ 10% Solution (Cat. No. 2127-50, 100, 500)

Octyl-β-D-glucopyranoside (Cat. No. 1646-5G, 25G, 50G)

Octyl-β-D-thioglucopyranoside (Cat. No. 1782-1G, 5G, 50G, 100G)

NP-40 (10% Solution in H<sub>2</sub>O) (Cat. No. 2111-100)

Pluronic® F-127 (Cat. No. 2730-10G, 50G)

Pluronic® F-127, MegaPure™ Detergent, 10% solution (Cat. No. 2731-25, 50, 100, 500))

SDS Solution (10% Solution in H<sub>2</sub>O) (Cat. No. 2102-100)

Sulfobetaine 3-16 (Cat. No. 1777-5G, 25G)

Taurocholic acid sodium salt (Cat. No. 2829-5G, 25G, 100G)

Triton X-100 (10% Solution in H<sub>2</sub>O) (Cat. No. 2104-100)

Triton X-100), MegaPure™ Detergent, 10% Solution (Cat. No. 2124-50, 100, 500)

Triton X-114), MegaPure™ Detergent, 10% Solution (Cat. No. 2137-50, 100, 500)

Tween® 20 (10% Solution in H<sub>2</sub>O) (Cat. No. 2111-100)

Tween® 20, MegaPure™, Detergent, 10% Solution (Cat. No. 2125-50, 100, 500)

Tween® 80, MegaPure™, Detergent, 10% Solution (Cat. No. 2126-50, 100, 500)

HANDLING: Do not take internally. Wear gloves and mask when handling

the product! Avoid contact by all modes of exposure.

FOR RESEARCH USE ONLY! Not to be used on humans.

BioVision Incorporated 155 S. Milpitas Boulevard, Milpitas, CA 95035 USA Tel: 408-493-1800 | Fax: 408-493-1801 www.biovision.com | tech@biovision.com