

Chemically Modified siRNA

(Cat# M1254-2 to M1254-10; 2 OD-10 OD; PAGE Purified; Store at -20°C)

I. Introduction:

BioVision offers low price and high quality **Chemically Modified siRNA Oligos**. All our siRNA oligos undergo vigorous process monitoring and strict quality control. Length and labeling are systematically controlled by PAGE or Mass Spectrometry Analysis. Quantity is systematically validated by UV absorbance at 260 nm. Although there are quite a lot principles for the designing of published sequences and quite a few companies supply on-line design service, up to now, no software can guarantee the efficiency of the oligos designed. BioVision and our collaboration laboratories have been involved in RNAi research for some time. BioVision's **chemically modified siRNA oligos** results in greater longevity in cell culture and stability in cell culture and in serum thereby enhancing the ability for this to be used for in vivo applications and has much more extensional effective time than the Standard siRNA.

II. Specifications of Chemically Modified siRNA Oligos:

- **Quality Control:** All our siRNA oligos undergo vigorous process monitoring and strict quality control. Length and labeling are systematically controlled by PAGE or Mass Spectrometry Analysis. Quantity is systematically validated by UV abs at 260 nm.
- **Purification:** Fully deprotected and desalted; **Purified by PAGE**
- **Length:** 19 to 23 mers
- **Bases:** RNA (A, C, G or U)
- **Backbone:** Phosphodiester bond
- **Labels and modifications:** Fluorescein, biotin and phosphate: 3' or 5' end
- **Format:** Single-strand RNA oligos is delivered in dry form
- **Oligonucleotide Technical Data Sheet:** Oligonucleotides are delivered with an Oligonucleotide Technical Data Sheet, which includes oligonucleotide name, sequence, concentration, size, purification method

III. Applications: siRNA Transfection

IV. Storage and Stability: Although oligonucleotides are stable in solution at 4°C for up to 2 weeks, we recommend storage should be at -20°C. Repetitive freeze-thaw cycles should be avoided by storing as aliquots. For long-term storage, siRNA oligos should be dried.

V. Shipment: Shipped by express delivery, dry in individual, transparent tubes

VI. Package Contents (Chemically Modified siRNA Oligo):

Components	Quantity	Part No.
Chemically Modified siRNA Oligo	Varies	M1254-XX-1
Universal Buffer	Varies	M1254-XX-2

VII. Related Products:

Product Name	Cat. No.	Quantity
Custom siRNA	M1253-2 to M1253-10	2 OD to 10 OD
Chemically Modified siRNA	M1254-2 to M1254-10	2 OD to 10 OD
Fluorescent labeled siRNA	M1255-2 to M1255-10	2 OD to 10 OD
Custom Neg Control siRNA	M1256-1	1 OD
Labeled Neg Control siRNA	M1257-1	1 OD
Custom Pos Control siRNA	M1258-1	1 OD
ss miRNA mimics	M1259-2 to M1259-10	2 OD to 10 OD
ss miRNA mimics Neg Control	M1260-1	1 OD
Labeled ss miRNA mimics Neg Control	M1261-1	1 OD
ds miRNA mimics	M1262-2 to M1262-10	2 OD to 10 OD
ds miRNA mimics Neg Control	M1263-1	1 OD
Labeled ds miRNA mimics Neg Control	M1264-1	1 OD
miRNA Inhibitor	M1265-2 to M1265-10	2 OD to 10 OD
miRNA Inhibitor Neg Control	M1266-1	1 OD
Labeled miRNA Inhibitor Neg Control	M1267-1	1 OD
Pre-designed siRNA Oligo Set A	M1268	Set
Pre-designed siRNA Oligo Set B	M1269	Set
Chemically Modified siRNA Set C	M1270	Set
Chemically Modified siRNA Set D	M1271	Set
Labeled siRNA Set E	M1272	Set
Labeled siRNA Set F	M1273	Set

FOR RESEARCH USE ONLY! Not to be used on humans.