

Wright-Giemsa Solution

CATALOG#:	M1438-500, -1000
AMOUNT:	-500 ml, -1000 ml
APPEARANCE:	Liquid
STORAGE:	Store at RT
DESCRIPTION:	The Wright-Giemsa Solution is intended for differential staining of blood smears, bone marrow and blood parasites (Cat# K1439). Use in conjunction with Phosphate Buffer Solution pH 6.8.
HANDLING:	For Laboratory Use Only. Not for drug, household or other uses. Do not take internally. Wear gloves and mask when handling the product! Avoid contact by all modes of exposure.

FOR RESEARCH USE ONLY! Not to be used in humans.