

980 Linda Vista Avenue Mountain View, CA 94043 USA Phone: (650)428-0236 Fax: (650)428-0336

PDG FRβ, Active

Recombinant protein expressed in Sf9 cells

Catalog # 7770-5 Lot#

Aliquot Size: 5 μg in 50 μl/vial

Shelf Life: 6-12 months from shipping date

Specific Activity: 20 nmol/min/mg

Product Description

Recombinant human PDGFR β (557-end) was expressed by baculovirus in Sf9 insect cells using a N-terminal GST tag. The gene accession number is NM 002609.

Gene Aliases

JTK12; PDGFR; CD140B; PDGFR1; PDGF-R-beta

Formulation

Recombinant protein stored in 50mM Tris-HCI, pH 7.5, 150mM NaCI, 0.25mM DTT, 0.1mM EGTA, 0.1mM EDTA, 0.1mM PMSF, 25% glycerol.

Storage and Stability

Store product at -70°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.

Scientific Background


PDGFRβ (platelet-derived growth factor receptor β) is a member of the PDGFR family of membrane receptors with intrinsic tyrosine kinase activity. PDGFRβ deficient mice are hemorrhagic, severely anemic and exhibit a defect in kidney glometric (1). However, absence of PDGFRβ has no impact on major blood vessels and the heart. PDGFRβ expression

Gentaur Europe BVBA Voortstraat 49, 1910 Kampenhout BELGIUM Tel 0032 16 58 90 45 info@gentaur.com

кетегепсеѕ

- Soriano, P: Abnormal kidney development and hematological disorders in PDGF beta-receptor mutant mice. Genes Dev. 1994 Aug 15;8(16):1888-96.
- 2. Xu, L et al: Blocking platelet-derived growth factor-D/platelet-derived growth factor receptor beta signaling inhibits human renal cell carcinoma progression in an orthotopic mouse model. Cancer Res. 2005 Jul 1;65(13):5711-9.

Specific Activity



The specific activity of PDGFRβ was determined to be 20 nmol/min/mg asper activity assay protocol.

Purity



Activity Assay Protocol

Reaction Components

Active Kinase #7770-5

Active PDGFR β (0.1 μ g/ μ l) diluted with Kinase Dilution Buffer and assayed as outlined in sample activity plot. (Note: these are suggested working dilutions and it is recommended that the researcher perform a serial dilution of Active PDGFR β for optimal results).

Kinase Dilution Buffer, pH 7.2

Kinase Assay Buffer II diluted at a 1:4 ratio (5X dilution) with 50 ng/ μl BSA solution.

Kinase Assay Buffer II, pH 7.2

Buffer components: 25mM MOPS, 12.5mM β -glycerolphosphate, 20mM MgC1 $_2$, 25mM MnC1 $_2$, 5mM EGTA, 2mM EDTA. Add 0.25mM DTTto Kinase Assay Buffer prior to use.

[32P]-ATP Assay Cocktail

Prepare 250 μ M [32 P]-ATP Assay Cocktail in a designated radioactive working area by adding the following components: 150 μ l of 10mM ATP Stock Solution, 100 μ l [32 P]-ATP (1mCi/100 μ l), 5.75ml of Kinase Assay Buffer. Store 1ml aliquots at -20° C.

10mM ATP Stock Solution

Prepare ATP stock solution by dissolving 55mg of ATP in 10ml of Kinase Assay Buffer. Store $200\mu l$ aliquots at $-20^{\circ}C$.

Substrate

Poly (Glu:Tyr, 4:1) synthetic peptide substrate diluted in distilled H_2O to a final concentration of 1 mg/ml.

Assay Protocol

Gentaur Europe BVBA Voortstraat 49, 1910 Kampenhout BELGIUM Tel 0032 16 58 90 45 info@gentaur.com



- Step 3. In a pre-cooled microfuge tube, add the following reaction components bringing the initial reaction volume up to 20ul:
 - Component 1. 10µl of diluted Active PDGFR\u00e3.
 - Component 2. 10µl of 1 mg/ml stock solution of substrate
- Step 4. Set up the blank control as outlined in step 3, excluding the addition of the substrate. Replace the substrate with an equal volume of distilled H₂O.

- Step 5. Initiate the reaction by the addition of 5μl [³²P]-ATP Assay Cocktail bringing the final volume up to 25μl and incubate the mixture in a water bath at 30°C for 15 minutes.
- Step 6. After the 15 minute incubation period, terminate the reaction by spotting 20µl of the reaction mixture onto individual pre-cut strips of phosphocellulose P81 paper.
- Step 7. Air dry the pre-cut P81 strip and sequentially wash in a 1% phosphoric acid solution (dilute 10ml of phosphoric acid and make a 1L solution with distilled H₂O) with constant gentle stirring. It is recommended that the strips be washed a total of 3 intervals for approximately 10 minutes each.
- Step 8. Count the radioactivity on the P81 paper in the presence of scintillation fluid in a scintillation counter.
- Step 9. Determine the corrected cpm by removing the blank control value (see Step 4) for each sample and calculate the kinase specific activity as outlined below.

Calculation of [P³²]-ATP Specific Activity (SA) (cpm/pmol)

Specific activity (SA) = cpm for 5μl [32P]-ATP / pmoles of ATP (in 5μl of a 250μM ATP stock solution, i.e., 1250 pmoles)

Kinase Specific Activity (SA) (pmol/min/μg or nmol/min/mg)

Corrected cpm from reaction / [(SA of 32 P-ATP in cpm/pmol)*(Reaction time in min)*(Enzyme amount in μg or mg)]*[(Reaction Volume) / (Spot Volume)]



Gentaur Europe BVBA Voortstraat 49, 1910 Kampenhout BELGIUM Tel 0032 16 58 90 45 <u>info@gentaur.com</u>