# RioVision

# Phospho-Bad (Ser112) Polyclonal Antibody

CATALOG #: 3269-100

AMOUNT: 100 µg

LOT#:

HOST: Rabbit

# FORMULATION:

100 µg (0.2 mg/ml) peptide affinity purified rabbit anti-Phospho-Bad (Ser112) polyclonal antibody in phosphate buffered saline (PBS), pH 7.2, containing 50% glycerol, 1% BSA, 0.02% thimerosal.

SPECIES REACTIVITY: Human, Mouse, Rat

#### STORAGE CONDITIONS:

Store at -20°C. For long term storage, aliquot and freeze at -70°C. Avoid repeated freeze/defrost cycles.

#### **DESCRIPTION:**

Bad is a member of the Bcl-2 family protein that selectively dimerize with Bcl-x<sub>L</sub> and Bcl-2 resulting in cell death. Survival factors such as IL-3 can inhibit the apoptotic activity of Bad by activating intracellular signaling pathways that result in the phosphorylation of Bad at Ser112 and Ser 136. Phosphorylation at these sites results in the binding of Bad to 14-3-3 proteins and the inhibition of Bad binding to Bcl-x<sub>L</sub> and Bcl-2. Recently, Akt has also been shown to promote cell survival via its ability to phosphorylate Bad at Ser136.

#### SPECIFICITY:

The antibody detects 23 kDa Bad only when phosphorylated at Ser112.

# APPLICATION AND USAGE:

The antibody can be used for Western blot analysis (1-2 µg/ml) and immunoprecipitation (5-20 µg/ml). However, the optimal conditions should be determined individually.

#### RELATED PRODUCTS:

rev.08/05

# Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

# Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

# Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

#### Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- **DNA Damage Quantification Kit**
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assav Kits

#### Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assav Kits
- Beta-Secretase Activity Assay Kit

# Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

# Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

# Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies

Gentaur Europe BVBA Voortstraat 49, 1910 Kampenhout BELGIUM Tel 0032 16 58 90 45 info@gentaur.com