

## Phospho-MAPKAPK-2 Polyclonal Antibody

**CATALOG #**: 3434-100

**AMOUNT**: 100 μg

LOT#:

HOST: Rabbit

**IMMUNOGEN:** Synthetic peptide corresponding to residues surrounding

Thr334 of human MAPKAPK-2

#### FORMULATION:

100 μg (0.5 mg/ml) antigen affinity purified rabbit polyclonal antibody in phosphate buffered saline (PBS), pH 7.2, containing 50% glycerol, 1% BSA, 0.02% thimerosal.

SPECIES REACTIVITY: Human, Mouse

#### STORAGE CONDITIONS:

Store at -20°C. For long-term storage, aliquot and freeze at -70°C. Avoid repeated freeze/defrost cycles.

#### BACKGROUND:

MAP kinase-activated protein kinase-2 (MAPKAPK-2), also known as p45 hsp27 kinase, is a 45-54 kDa serine-threonine protein kinase. MAPKAPK-2 is rapidly phosphorylated and activated in response to stress, cytokines, and chemotactic factors. One of the major substrates of MAPKAPK-2 is hsp27, which stimulates actin polymerization in order to facilitate recovery from destruction of cytoskeleton during cellular stresses. Two isoforms are produced due to alternative splicing of the same gene which differs in their C-terminals. This antibody recognizes both isoforms.

#### APPLICATION AND USAGE:

The antibody can be used for Western blotting (1  $\mu$ g/ml). However, the optimal conditions should be determined individually.

#### SPECIFICITY:

The antibody detects phosphorylated MAPKAPK-2 (Thr334) of human and mouse origins. Reactivity to other species has not been tested.

# FOR RESEARCH LISE ONLY! Not to be used on humans

#### **RELATED PRODUCTS:**

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

### Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

#### Cell Proliferation & Senescence

- Quick Cell Proliferation Assav Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

## Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

#### Signal Transduction

- cAMP & cGMP Assav Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

#### Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

## Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

### **Growth Factors and Cytokines**

Monoclonal and Polyclonal Antibodies

