

# Phospho-MEK1/2 Polyclonal Antibody

CATALOG #:	3519-100
AMOUNT:	100 μΙ
_OT #:	
HOST:	Rabbit

#### FORMULATION:

100 µl purified rabbit polyclonal antibody in phosphate-buffered saline (PBS) containing 50% glycerol, 1% BSA, and 0.02% sodium azide.

SPECIES REACTIVITY: Human, Mouse, Rat

#### STORAGE CONDITIONS:

Store at -20°C. For long term storage, aliquot and freeze at -70°C. Avoid repeated freeze/defrost cycles.

#### BACKGROUND:

MAP kinase kinase 1 (Mek1), also known as mitogen activated protein kinase kinase, MAP kinase kinase, ERK kinase 1 or microtubule-associated protein 2 kinase, is a 43 kDa dual-specificity protein kinase that phosphorylates threonine and tyrosine residues on MAP kinases such as Erk 1 and 2. Mek1 and Mek2 are involved in a diverse array of cellular processes such as stress-activated response, apoptosis, cytokine-induced cell proliferation, and DNA recombination during meiosis. Mek1 and Mek2 are activated by phosphorylation on serine 218 and 222 by serine-threonine kinase, Raf-1, which is part of the p21ras signal transduction pathway.

### SPECIFICITY:

The antibody recognizes 45 kDa phospho-Mek1/2 (Ser217/221). It does not detect the unphosphorylated protein.

#### APPLICATION:

The antibody can be used in Western blotting (1:200-1000) and immunocytochemistry (1:100-200). However, the optimal concentrations should be determined individually.

FOR RESEARCH USE ONLY! Not to be used on humans.

#### RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

## Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

#### Cell Proliferation & Senescence

- · Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

## Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assavs & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

### Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

## Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

## Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies

