

Viral Transport, Veal Infusion Medium

04/21

CATALOG #:

M1516-50	50 ml
M1516-100	100 ml
M1516-200	2 x 100 ml
M1516-500	500 ml

FORM: Clear liquid with a light amber tint

FEATURES:

- Ready-to-use sterile solution
- For molecular assays and rapid antigen testing
- Direct fluorescent antibody applications

FORMULATION: In solution containing a veal infusion broth, Bovine Serum Albumin, Gentamicin and Amphotericin.

DESCRIPTION: **BioVision's Viral Transport, Veal Infusion Medium** is formulated for viral research and diagnostics. It is used for the collection, transport and culture of virus. Its unique formulation is negative for microbial growth after 48 h at 37 °C. This medium can be used for multiple applications such as antigen testing, fluorescence assays, etc.

STORAGE: Store at 4 °C. The product is stable for one year when stored as recommended.

HANDLING: Do not take internally. Wear gloves and mask when handling the product! Avoid contact by all modes of exposure.

RELATED PRODUCTS:

BioVision Product Name	Cat. No.	Size
Viral Transport Medium	M1515	50, 100, 2 x 100, 500 ml
Mag-Lentivirus and Retrovirus Purification Kit	K1458	20, 100 Preps
Adenovirus Mini Purification Kit	K1300	10, 20 Preps
Lentivirus Mini Purification Kit	K1305	10, 20 Preps
Retrovirus Mini Purification Kit	K1307	10, 20 Preps
Viral DNA extraction kit	K1446	100 Preps
PEG Virus Precipitation Kit	K904	50, 200 preps
Coronavirus Enhanced Rapid RT-qPCR Detection Kit	K1484	100 Rxns

DISCLAIMER: *FOR RESEARCH USE ONLY! Not to be used on humans.*