

Eotaxin/CCL11, Human Recombinant

CATALOG #:	4028-20	20 µg
	4028-100	100 µg
	4028-1000	1 mg
SOURCE:	<i>E. coli</i>	
PURITY:	>97% by SDS-PAGE and HPLC analyses Endotoxin level is <0.1 ng per µg of Eotaxin	
MOL. WEIGHT:	8.3 kDa	
FORM:	Lyophilized from 20 mM phosphate buffer (pH 7.4) and 150 mM NaCl.	
SPECIFIC ACTIVITY:	>1x10 ⁵ IU/mg.	

RECONSTITUTION:

The lyophilized Eotaxin can be reconstituted in water to a concentration of 0.1-1 µg/µl. It is recommended that further dilutions be made into buffer containing carrier protein or medium containing serum.

STORAGE CONDITIONS:

The lyophilized Eotaxin is best-stored desiccated below 0°C. Reconstituted Eotaxin should be stored at 4°C or in working aliquots at -20°C. For long term storage add 0.1% BSA as a carrier protein

DESCRIPTION:

Eotaxin is a recently discovered chemokine characterized by its high chemotactic selectivity for eosinophils. The position of the four cystein residues in Eotaxin places it in the C-C family of chemokines, along with chemokines such as RANTES, MCP-3, and MIP-1 alpha. The human Eotaxin is an 8.3 kDa protein containing 74 amino acid residues. The sequence of the first five N-terminal amino acids is Gly-Pro-Ala-Ser-Val.

BIOLOGICAL ACTIVITY:

The ED₅₀ of human Eotaxin by chemotaxis based assays is 10 ng/ml using human peripheral blood eosinophils.

FOR RESEARCH USE ONLY! Not to be used in humans

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines Monoclonal and Polyclonal Antibodies