BioVision

Eotaxin/CCL11, Murine Recombinant

CATALOG #: 4029-10 10 μg

4029-1000 1 mg

LOT #: _____

SOURCE: E. coli

MOL. WEIGHT: 8.4 kDa

PURITY:

>98% by SDS-PAGE and HPLC analyses. Under reducing SDS-PAGE conditions the apparent molecular weight of murine eotaxin is as expected (8.4 kDa), while under nonreducing conditions the apparent molecular weight is 17 kDa. Endotoxin level is <0.1 ng per μ g of murine eotaxin.

FORM:

Sterile filtered and then lyophilized with no additives. The protein may appear as a haze or film that is difficult to see at the bottom of the vial.

RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute in water to a concentration of 0.1-1 μ g/ μ l. It is recommended that further dilutions be made into buffer containing carrier protein or medium containing serum.

STORAGE CONDITIONS:

The lyophilized eotaxin is best-stored desiccated below 0°C. Reconstituted eotaxin should be stored at 4°C or in working aliquots at -20°C.

DESCRIPTION:

Murine eotaxin shares many of the biological and physical characteristics of human eotaxin. The chemokine exhibits a high chemotactic selectivity for eosinophils and belongs to the C-C family of chemokines, along with chemokines such as RANTES, MCP-3, and MIP-1 alpha. Murine eotaxin is an 8.4 kDa protein containing 74 amino acid residues.

BIOLOGICAL ACTIVITY:

The recombinant murine eotaxin was found to induce chemotaxis of purified eosinophils at concentrations ranging between 100 to 1000 ng/ml.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies

