

Recombinant Human Epiregulin (EREG)

CATALOG #: 4069-25 25 µg
4069-1000 1 mg

LOT #: _____

SOURCE: *E. coli*

PURITY: >97% by SDS-PAGE and HPLC analyses
Endotoxin level is <0.1 ng per µg of EREG.

MOL. WEIGHT: 6.0 kDa

FORM: Sterile filtered and lyophilized from 0.5 mg/ml solution containing 20 mM PBS, pH 7.4 and 130 mM NaCl.

RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute in sterile water to 0.1–1 mg /ml. The solution can then further diluted to other aqueous solutions.

STORAGE CONDITIONS:

The lyophilized EREG is best-stored desiccated below 0°C. Reconstituted EREG should be stored at working aliquots at –20°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).

DESCRIPTION:

Epiregulin is part of the EGF family. Epiregulin functions as a ligand of EGFR, as well as a ligand of most members of the ERBB (v-erb-b2 oncogene homolog) family of tyrosine-kinase receptors. Epiregulin is expressed mostly in the placenta and peripheral blood leukocytes and in specific carcinomas of the bladder, lung, kidney and colon. Epiregulin stimulates the proliferation of keratinocytes, hepatocytes, fibroblasts and vascular smooth muscle cells. Epiregulin inhibits the growth of several tumor-derived epithelial cell lines. Human Epiregulin is initially synthesized as a glycosylated 19.0 kDa transmembrane precursor protein, which is processed by proteolytic cleavage to produce a 6.0 kDa mature secreted sequence. Recombinant human EREG produced in *E.Coli* is a single, non-glycosylated, polypeptide chain containing 50 amino acids and having a molecular mass of 6 kDa.

BIOLOGICAL ACTIVITY:

The ED₅₀ as determined by the dose-dependent stimulation of the proliferation of murine balb/3T3 cells is ≤2.0 ng/ml, corresponding to a specific activity of ≥5 x 10⁵ units/mg.

FOR RESEARCH USE ONLY! Not to be used on humans.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies