# BioVision

#### rev. 07/07

# **Recombinant Murine MCP-3**

CATALOG #:	4229-10 4229-1000	10 µg 1 mg
LOT #:		
SOURCE:	E. coli	
PURITY:	>99% by SDS-PAGE and HPLC analyses Endotoxin level is <0.1 ng per µg of MCP-3.	
MOL. WEIGHT:	8.5 kDa	
FORM:	Sterile filtered and lyophilized with no additives	

## **RECONSTITUTION:**

Centrifuge the vial prior to opening. Reconstitute in water to a concentration of 0.1-1 mg/ml. This solution can then be diluted into other aqueous buffers or stored at  $4^{\circ}$ C for 1 week or – 20°C for future use.

#### STORAGE CONDITIONS:

The lyophilized murine MCP-3 is best-stored desiccated below 0°C. Reconstituted MCP-3 should be stored in working aliquots at -20°C.

#### **DESCRIPTION:**

Murine MCP-3 (Macrophage/Monocyte Chemotactic Protein-3) plays an important role in the inflammatory response of blood monocytes and T-cells. The recombinant murine MCP-3 is a 8.5 kDa protein containing 74 amino acid residues.

## **BIOLOGICAL ACTIVITY:**

The biological activity of murine MCP-3 was determined by a chemotaxis assay using Balb/C mouse spleen MNCs using a concentration range of 10-100 ng/ml.

# RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits
  Signal Transduction

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit
- Molecular Biology & Reporter Assays
  - siRNA Vectors
  - Cloning Insert Quick Screening Kit
  - Mitochondrial & Genomic DNA Isolation Kits
  - 5 Minutes DNA Ligation Kit
  - 20 Minutes Gel Staining/Destaining Kit
  - β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines Monoclonal and Polyclonal Antibodies



Gentaur Europe BVBA Voortstraat 49, 1910 Kampenhout BELGIUM Tel 0032 16 58 90 45 info@gentaur.com

FOR REAFARAULUAE AND VENE ( - L. ...........