

# **Recombinant Murine MIP-2**

**CATALOG #**: 4251-20 20 μg

4251-1000 1 mg

LOT #: \_\_\_\_\_

SOURCE: E. coli

**PURITY:** >98% by SDS-PAGE and HPLC analyses

Endotoxin level is <0.1 ng per μg of MIP-2.

MOL. WEIGHT: 7.8 kDa

**FORM:** Sterile filtered and then lyophilized with no additives.

### RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute in water to a concentration of 0.1-1 mg/ml. The solution can be diluted into other aqueous buffers and store at 4°C for 1 week or –20°C for future use.

# STORAGE CONDITIONS:

The lyophilized MIP-2 is best-stored desiccated below 0°C. Reconstituted murine MIP-2 should be stored in working aliquots at -20°C. Avoid freeze/thaw cycles.

#### DESCRIPTION:

Murine microphage inflammatory protein-2 (MIP-2) promotes neutrophil chemotaxis and degranulation. BioVision's murine MIP-2 is a 7.8 kDa protein containing 73 amino acid residues including the 'RLR' motif common to the CXC chemokine family that bind to CXCR1 or CXCR2.

### **BIOLOGICAL ACTIVITY:**

Determined by measuring the dose dependent mobilization of intracellular calcium (calcium flux) with human neutrophils. Significant calcium mobilization is observed with  $\geq 5$  ng/ml of recombinant mouse MIP-2.

# **RELATED PRODUCTS:**

# Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

# Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

### Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

### Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

### Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

# Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

# Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

# **Growth Factors and Cytokines**

Monoclonal and Polyclonal Antibodies

