# **BioVision**

# **Recombinant Human Midkine**

**CATALOG #**: 4479-20 20 μg

4479-100 100 μg 4479-1000 1 mg

LOT #: \_\_\_\_\_

SOURCE: E. coli

**PURITY:** >98% by SDS-PAGE and HPLC analyses

Endotoxin level is <0.1 ng per μg.

MOL. WEIGHT: 13.4 kDa

**FORM:** Lyophilized from 0.4X PBS, pH7.4

#### RECONSTITUTION:

Centrifuge the vial before opening. Reconstitute in water to a concentration of 0.1-1.0 mg/ml. This solution can then be diluted into other aqueous buffers and stored at 4°C for 1 week or – 20°C for future use.

#### STORAGE CONDITIONS:

The lyophilized protein is best-stored desiccated below 0°C. Reconstituted human Midkine should be stored in working aliquots at -20°C.

# **DESCRIPTION:**

Human Midkine (MK) is a new member of the heparin-binding neurotrophic factor family. MK and PTN (Pleiotrophin) are structural homologs, and are highly conserved among species. MK plays important roles in development and carcinogenesis and has several important biological effects, including promotion of neurite extension and neuronal survival. Recombinant midkine is a 13.4 kDa protein, comprising of 123 amino acid residues.

## **BIOLOGICAL ACTIVITY:**

Determined by its ability to chemoattract human neutrophils using a concentration of 0.1-10 ng/ml.

#### RELATED PRODUCTS:

## Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

#### Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

#### Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

# Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

# Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

#### Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

#### Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

#### **Growth Factors and Cytokines**

