BioVision

Recombinant Human Eotaxin-2

CATALOG #: 4497-20 20 μg

4497-1000 1 mg

LOT #: _____

SOURCE: E. coli

PURITY: >97% by SDS-PAGE analyses

Endotoxin level is <0.1 ng per µg.

MOL. WEIGHT: 8.8 kDa

FORM: The protein was lyophilized from a concentrated (1 mg/ml)

sterile solution containing 20 mM PBS pH-7.4 and 0.15 M

sodium chloride.

RECONSTITUTION: Centrifuge the vial prior to opening. Reconstitute in water to

a concentration of 0.1-1.0 mg/ml. This solution can be stored at 4°C for 1 week or –20°C for future use. It is recommended that further dilutions be made into buffer containing carrier

protein or medium containing serum.

STORAGE CONDITIONS: The lyophilized protein is best-stored desiccated below 0°C.

Reconstituted human Eotaxin-2 should be stored in working aliquots at -20°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent

freeze-thaw cycles.

DESCRIPTION: Eotaxin-2 also known as MPIF-2 is a CC chemokine

characterized by its high chemotactic selectivity for eosinophils. BioVision's human eotaxin-2 is 8.8 kDa protein

containing 78 amino acid residues.

BIOLOGICAL ACTIVITY: Determined by its ability to chemoattract using peripheral

blood eosinophils using a concentration of 50-100 ng/ml.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

