# BioVision

# Recombinant Murine IL-22

CATALOG #:	4602-10 4602-50 4602-1000	10 μg 50 μg 1 mg
LOT #:		
SOURCE:	E. coli	
PURITY:	>97% by SDS-PAGE and HPLC Endotoxin level is <0.1 ng per µg of IL-22.	
MOL. WEIGHT:	33.4 kDa	
FORM:	Lyophilized from 1X PBS Solution	

#### **RECONSTITUTION:**

Centrifuge the vial prior to opening. Reconstitute in water to a concentration of 0.1 to 1.0 mg/ml. This solution can then be stored at 4°C for 1 week or -20°C for future use. For long term storage, it is recommended to add a carrier protein (e.g., 0.1% BSA). Avoid repeating freeze/thaw cycles.

# STORAGE CONDITIONS:

The lyophilized IL-22 is best-stored desiccated below 0°C. Reconstituted IL-22 should be stored at working aliquots at -20°.

#### **DESCRIPTION:**

Murine Interleukin-22 (IL-22), also known as T cell-derived inducible factor (ILTIF), is a cytokine that is structurally related to IL-10. Studies have shown that IL-22 induces acute-phase reactants *in vitro* and *in vivo*, suggesting that the cytokine contributes to the regulation of inflammatory responses. Recombinant murine IL-22 contains 2 x 146 amino acid chains having a total molecular mass of 33.4 kDa.

# ACTIVITY:

Recombinant murine IL-22 is fully biologically active when compared to a standard. The ED50 as determined by its ability to induce IL-10 secretion in Colo205 cells is less than 0.5 ng/ml, corresponding to a specific activity of  $2 \times 10^6$  IU/mg.

#### **RELATED PRODUCTS:**

rev. 03/09

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit

Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits
Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit
- Molecular Biology & Reporter Assays
  - siRNA Vectors
  - Cloning Insert Quick Screening Kit
  - Mitochondrial & Genomic DNA Isolation Kits
  - 5 Minutes DNA Ligation Kit
  - 20 Minutes Gel Staining/Destaining Kit
  - β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

# Growth Factors and Cytokines

Manageral and Debulanet Astheadies

