

# **Recombinant Mouse Growth Hormone (GH)**

**CATALOG #**: 4770-50 50 μg

4770-1000 1 mg

LOT #: \_\_\_\_\_

SOURCE: E. coli

**PURITY:** ≥ 95% as determined by SDS-PAGE and HPLC

**ENDOTOXIN:** Less than 0.1 ng/mg of GH

PHYSICAL APPEARANCE: Lyophilized powder

# **RECONSTITUTION:**

Centrifuge the vial prior to opening. It is recommended to reconstitute the lyophilized mGH in sterile deionized water to 0.1-1.0 mg/ml. This solution can then be further diluted to other aqueous solutions or store at -20C for future use.

# STORAGE CONDITIONS:

Lyophilized GH should be stored desiccated below 0°C. Upon reconstitution GH can be stored at 4°C for 2-7 days. For long-term storage, it is recommended to add a carrier protein (0.1% HSA or BSA) and store aliquots at –20°C or –70°C. Avoid freeze-thaw cycles.

### DESCRIPTION:

GH is a member of the somatotropin/prolactin family of hormone which plays an important role in growth control. Recombinant mouse Growth Hormone produced in *E.Coli* is a single, non-glycosylated, single chain containing 191 amino acids and having a molecular mass of 22.0 kDa.

### **BIOLOGICAL ACTIVITY:**

Recombinant mouse Growth Hormone is fully biologically active when compared to standard human growth hormone which is 3 units/mg.

# **RELATED PRODUCTS:**

### Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

# Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

### Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

### Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

# Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

# Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

# Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

County Forton and Outstand Management and Delivational Authorities

Gentaur

EOR RESEARCH LISE ONLY I Not to be used on humans