

## Human recombinant linear tetra-ubiquitin


<b>CATALOG #:</b>	6426-100	100 µg
<b>ALTERNATE NAMES:</b>	UB4	
<b>MOL. WEIGHT:</b>	34.8 kDa (His Tagged)	
<b>FORMULATION:</b>	In 20 mM Tris-HCl, pH 7.5, 0.15 M NaCl and 1 mM EDTA	
<b>STORAGE CONDITIONS:</b>	Aliquot and store at -80°C. Avoid repeated freezing and thawing cycles.	

### DESCRIPTION:

A wide range of cellular processes are modulated through the generation and attachment of polyubiquitin (polyUb) chains to target proteins. Increasing evidence suggests that polyUb chains joined through linear peptide bonds between the C-terminus of a ubiquitin and the N-terminus of another play important functional roles. The enzyme machinery responsible for the generation of linear polyUb chains has been termed LUBAC, consisting of HOIL-1L and HOIP. Chains of these type have been determined to have an open conformation, similar to polyUb K63, but with very distinct functional properties. Linear polyUb chains are cleaved by the deubiquitylases CYLD, USP5 (IsoT), USP2 and have been shown to bind to many UBDs including NEMO and Trabin-n (3xznf). Recombinant linear chains of defined length are expressed in *E. coli* and purified to homogeneity. Amide linkages join the N- and C-terminus of each ubiquitin molecule to each other. This molecule is HIS-tagged at the N-terminus of the most distal ubiquitin.

### APPLICATIONS:

- Investigation of binding interactions through pull-down studies
- Determining potential deubiquitylase activity towards linear chain



Approximately 50ng of polyUb chains were subjected to SDS-PAGE analysis (15%) and visualized by silver staining.

### RELATED PRODUCTS:

- K48-linked di-ubiquitin (Cat. No. 6415-50)
- K48-linked tri-ubiquitin (Cat. No. 6416-50)
- K48-linked tetra-ubiquitin (Cat. No. 6417-25)
- K63-linked di-ubiquitin (Cat. No. 6418-50)
- K63-linked tri-ubiquitin (Cat. No. 6419-50)
- K63-linked tetra-ubiquitin (Cat. No. 6420-25)
- K11-linked di-ubiquitin (Cat. No. 6421-50)
- K11-linked tri-ubiquitin (Cat. No. 6422-50)
- K11-linked tetra-ubiquitin (Cat. No. 6423-25)
- Human recombinant Linear di-ubiquitin (Cat. No. 6424-100)
- Human recombinant Linear tri-ubiquitin (Cat. No. 6425-100)
- Human recombinant Linear penta-ubiquitin (Cat. No. 6427-100)
- Di-ubiquitin explorer panel (Cat. No. 6428-5)

**FOR RESEARCH USE ONLY! Not to be used in humans.**