

Recombinant c-Jun Fusion Protein

CATALOG #: 7001-100

AMOUNT: 100 µg (2 mg/ml)

LOT #: _____

BACKGROUND:

c-Jun is a component of the transcription factor AP-1 that binds and activates transcription at TRE/AP-1 elements. The transcriptional activity of c-Jun is regulated by phosphorylation at Ser63 and Ser73. Extracellular signals including growth factors, transforming oncoproteins, and UV irradiation stimulate phosphorylation of c-Jun at Ser63/73 and activate c-Jun dependent transcription. The MAP kinase homologue, JNK, binds to the N-terminal region of c-Jun and phosphorylates c-Jun at Ser63/73. The activity of JNK is stimulated by the same signals that activates c-Jun.

PREPARATION:

c-Jun fusion protein was expressed in *E. coli*. It is expressed as a recombinant protein fusion of human c-Jun residue 1-79 and glutathione S-transferase. The protein was then purified by affinity chromatography.

PURITY: ≥ 90%

FORM:

Liquid. 2 mg/ml in storage buffer containing 20 mM Tris-HCl (pH 7.5), 50 mM NaCl, 2 mM EDTA, 1 mM DTT, and 50% glycerol.

MOL. WEIGHT: 37 kDa

STORAGE: -20°C

ACTIVITY:

C-Jun fusion protein at a concentration of 2 µg/20 µl reaction was phosphorylated in an *in vitro* kinase assay with 200 µM ATP. After 30 minutes incubation at 30°C, phosphorylation was detected by Western blot with phospho-specific c-Jun antibodies.

REFERENCES:

Binetruy, *et al.* (1991) *Nature* 351:122-127.
Smeal, T. *et al.* (1991) *Nature* 354, 494-496.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay Kit
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Quality Antibodies for Apoptosis and Signal Transduction Molecules